



# A Smarter Solution

## Syntri Safeguard

The accuracy of patient tests often hinges upon proper preparation and consistency in processes throughout the laboratory. When paired with Syntri Service and reagents featuring Syntri Safeguard RFID technology, the Revos tissue processor tracks and logs reagent usage. It remembers every bottle loaded, protocol executed and service performed – helping to reduce downtime and facilitating efficient use of reagents.



Find it at [fisherhealthcare.com](http://fisherhealthcare.com)

 **fisher**  
healthcare  
part of Thermo Fisher Scientific



For more information about the Revos tissue processor, please visit [epredia.com/revos](https://epredia.com/revos)

### Syntri Safeguard Consumables

Description	Cat. Number
Syntri Safeguard™ 10% Neutral Buffered Formalin, 4 x 5 L	5735SS
Syntri Safeguard™ Pen-Fix™, 4 x 5 L	6115SS
Syntri Safeguard Dehydrant 100%, 4 x 5 L	6215SS
Syntri Safeguard Xylene, 4 x 5 L	6615SS
Syntri Safeguard™ Clear-Rite™ 3, 4 x 5 L	6915SS
Syntri Safeguard™ Flex™ 100, 4 x 5 L	8115SS
Syntri Safeguard Isopropyl Alcohol, 4 x 5 L	9517SS
Syntri Safeguard™ Formal-Fixx™, 4 x 5 L	6764215SS
Shandon™ Syntri Safeguard™ Xylene Substitute, 4 x 5 L	6764515SS

Distributed by Fisher Healthcare

Contact us today:

#### In the United States

Order online: [fisherhealthcare.com](https://fisherhealthcare.com)

Fax an order: 1-800-290-0290

Call customer service: 1-800-640-0640



© 2021 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at [fishersci.com/trademarks](https://fishersci.com/trademarks).

BN20215914