

Nothing has changed, except everything

Introducing the SeqStudio Genetic Analyzer

Same workflow, same trusted technology—now with an innovative all-in-one cartridge that reduces setup time from hours to minutes. The new Applied Biosystems™ SeqStudio™ Genetic Analyzer is a low-throughput, easy-to-use, and convenient benchtop system that delivers gold-standard Sanger sequencing and fragment analysis with just a simple click. It is easily used across a broad range of applications as well as by multiple users.

Just click. That's it.

Load up your samples, click in the cartridge, and go—no need for technical experts for setup or maintenance, no calibration before each run, and no installation of capillaries, polymer, or anode buffer. The cartridge is stable for four months on the instrument.

Access your data anytime, anywhere*

This cloud-enabled system is integrated with Thermo Fisher Connect so you can set up your run, check its progress, and access your data from wherever you have Internet access, all without the need to stay in the lab.* Cloud-based Applied Biosystems™ Sanger analysis modules are available for secondary analysis.

Get the same data quality, service, and support you've come to expect from Applied Biosystems™ technology with a modernized experience at a price any lab can afford.



- 

Convenient and easy-to-use cartridge system includes the capillaries, polymer, and buffers required for each run
- 

Combine Sanger sequencing and fragment analysis reactions in the same instrument run
- 

Fast turnaround with a run time of as little as 30 minutes
- 

Set up and monitor your run. View, manage, and share your data from anywhere, using your phone, tablet, or computer (Mac™ computer or PC).*
- 

Secondary analysis software package included with system purchase
- 

One-year instrument warranty included
- 

Get up and running quickly with SmartStart orientation—in-lab, one-day training for basic setup, cloud enablement, reagent review, and instrument operation and maintenance

* Internet connection and Thermo Fisher Cloud account required.

Same workflow, same trusted data. Now with one-click convenience.

Sanger sequencing workflow

Primer design and amplification



Amplify the DNA region of interest: select from pre-designed primer pairs or design your own

PCR cleanup



Clean up your sample with one-tube, one-step PCR clean-up using Applied Biosystems™ ExoSAP-IT™ Express reagent

Cycle sequencing



Generate sequencing products by cycle sequencing

Sequencing cleanup



Purify reactions to remove unincorporated Applied Biosystems™ BigDye™ terminators and salts

Analyze Sanger sequencing data



Translate raw data into corresponding electropherograms

Fragment analysis workflow

PCR amplification



Amplify the DNA fragment of interest by PCR using 5' fluorescent-labeled primer

Prepare sample for analysis



Add aliquot of PCR product to size standard, denature, and proceed with electrophoresis

Run samples

Perform capillary electrophoresis

Analyze fragment analysis data



Analyze the fragment profile

Applications

Sanger sequencing

Resequencing, NGS confirmation, indels, heterozygote detection, and minor variant detection

Fragment analysis

Microsatellite analysis, compatible with MLPA™ amplification, Applied Biosystems™ SNaPshot™ applications, and cell line authentication

Service, support, and training

We offer a range of service and training options to help maintain instrument performance, helping you achieve your research outcomes. Every new SeqStudio instrument comes with SmartStart support to get you up and running quickly and covers basic setup, cloud enablement, printer networking, starter reagent review, software use, and instrument operation and maintenance.

Ordering information

Product	Cat. No.
SeqStudio Genetic Analyzer System with SmartStart Includes: <ul style="list-style-type: none"> SeqStudio Genetic Analyzer SeqStudio Genetic Analysis Software SmartStart SeqStudio Full Day 1-year warranty 	A35644
SeqStudio Genetic Analyzer System with SmartStart plus 1-year extended warranty <ul style="list-style-type: none"> Includes all items from A35644 plus additional 1-year warranty 	A35645
SeqStudio Genetic Analyzer System with SmartStart plus 3-year extended warranty <ul style="list-style-type: none"> Includes all items from A35644 plus additional 3-year warranty 	A35646
SeqStudio Starter Kit	A35000
SeqStudio Cartridge	A33671
Cathode Buffer Container (4 pk)	A33401
Cathode Buffer Container Reservoir Septa (20 pk)	A35640
96-Well Plate Septa (20 pk)	A35641
8-Strip Tube Septa (24 pk)	A35643
Integrated Capillary Protector	A31923

Find out more at thermofisher.com/seqstudio

In the United States:

For customer service, call 1-800-766-7000

To fax an order, use 1-800-926-1166

To order online: fishersci.com

In Canada:

For customer service, call 1-800-234-7437

To fax an order, use 1-800-463-2996

To order online: fishersci.ca



part of Thermo Fisher Scientific

For Research Use Only. Not for use in diagnostic procedures. © 2018 Thermo Fisher Scientific Inc. All rights reserved.

All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Mac is a trademark of Apple Inc. Multiplex Ligation-dependent Probe Amplification (MLPA) is a trademark of MRC-Holland. **BN20170956**

COL15038 1017