



Because every section matters

Epredia Ultra Microtome Blades

Obtaining high-quality sections and minimizing wasted tissue relies on high-quality microtomes paired with high-quality microtome blades. Epredia™ Ultra microtome blades offer high consistency, impressive durability, dependability, and precise performance.

Our proprietary manufacturing process employs a revolutionary blade coating system shown to increase strength and durability of microtome blades. Our proprietary steel hardening technology (Pink Technology) offers outstanding cutting performance and longevity.

The Ultra Advantage

- High- and low-profile blades are compatible with most manufacturers' microtomes and cryostats
- Angle profile facilitates fast ribboning
- Proprietary manufacturing process produces an edge twice as strong as traditional blades
- Special coating treatment provides durability and consistency from blade to blade



Product	Cat. No.
Epredia MX35 Ultra low-profile blades (50 per pack)	3053835
Epredia HP35 Ultra high-profile blades (50 per pack)	3153735

Distributed by Fisher Healthcare

Contact us today:

In the United States

Order online: fisherhealthcare.com

Fax an order: 1-800-290-0290

Call customer service: 1-800-640-0640

