

DWK LIFE SCIENCES INTRODUCES **KIMBLE® SILICONE LIDS** SPLASH YOUR LAB WITH **COLOR!**

**NEW
PRODUCT**

★★★★★

The reusable KIMBLE® silicone lids stretch easily and adapt to many container sizes, improves lab safety, enhances sample management, offers a more sustainable and cost-effective solution over foils, films and other single-use covers.

VERSATILE

- Fits on beakers, flasks, bottles, cylinders, plastic and glass products
- Each size conforms to various diameters
- Works with spouts, beads and thread finishes
- Adaptable for many uses beyond limits of traditional films and covers

LAB SAFETY

- Improves safety performance over foils or films
- Color aids quick visual identification
- Protects skin and nose from contents
- Defends against spill/leakage

SAMPLE MANAGEMENT

- Reduces cross-contamination of samples
- Compatible with lab markers to enhance sample identification
- Improves your lab workflow efficiency
- Uses color to track projects and samples

SUSTAINABLE

- Reusable, customizable and dishwasher safe
- Reduces non-reusable consumables that generate waste such as films and foils
- Chemically resistant, autoclavable and pliable over a wide temperature range
- Reduces environmental impact and lowers lab operational costs



KIMBLE® SILICONE LID

Three flexible sizes, three bright colors, versatile and simple to use: The new way to cover vessels.

HANDLING THE SILICONE LID



Position the edge of the vessel opening in the groove located on the inner surface of the lid.



Now hold the silicone lid onto the outer edge with one hand and pull the rest of the lid with the other hand over the entire opening.



To ensure that the silicone lid fits correctly, it must be tightly seated on the vessel.

MATERIAL RESISTANCE

| Substance groups + 23 °C | Silicone |
|---------------------------|----------|
| Acetone | ++ |
| Acetonitrile | + |
| Chloroform | ++ |
| Dichloromethane | ++ |
| Dimethyl formamide (DMF) | + |
| Dimethyl sulfoxide (DMSO) | ++ |
| Ether (Diethylether) | ++ |
| Ethyl Alcohol | ++ |
| Hexane | ++ |
| Isopropyl Alcohol | ++ |
| Methyl Alcohol | ++ |
| Tetrahydrofuran (THF) | ++ |
| Toluene | + |

++ = very good resistance

+ = good to conditional resistance

| Size | Opening diameter range | Qty Pack | Catalog no. | | |
|------------------|--------------------------|----------|-------------|------------|------------|
| | | | Pink | Cyan | Green |
| S | Stretch Ø 43 – 61 mm | 1 | 15-112-171 | 15-112-172 | 15-112-173 |
| M | Stretch Ø 64 – 76 mm | 1 | 15-112-174 | 15-112-175 | 15-112-176 |
| L | Stretch Ø 84 – 116 mm | 1 | 15-112-177 | 15-112-178 | 15-112-179 |
| Set S-M-L | 1 piece each size | 3 | 15-112-180 | 15-112-181 | 15-112-182 |

TECHNICAL PRODUCT INFORMATION

The KIMBLE® silicone lid is made of stretchable, chemically resistant and heat-resistant silicone. The recommended usable temperature range is from -40 to +180 °C.

Note: Please note that prolonged exposure to a solvent will cause swelling of the material. Product applications with direct solvent contact must be tested and evaluated by the user prior to the start of the test. Furthermore, the relevant health and safety regulations must be observed.

KIMBLE® SILICONE LIDS

GLASSWARE COMPATIBILITY

BEAKERS

| Cat. No. | Description | Capacity (mL) | Size Compatibility | Case Qty |
|-----------------------------------------|------------------|---------------|--------------------|----------|
| KIMBLE® KIMAX® Colorware Beakers | | | | |
| 10-200 | Blue | 50 | Small | 12 |
| 10-200-50 | Green | 50 | Small | 12 |
| 10-200-51 | Red | 50 | Small | 12 |
| 02-555-107 | Pink | 50 | Small | 12 |
| 10-200-3 | Yellow | 50 | Small | 12 |
| 10-200-48 | Blue | 100 | Small | 12 |
| 10-200-49 | Green | 100 | Small | 12 |
| 10-200-6 | Red | 100 | Small | 12 |
| 02-555-108 | Pink | 100 | Small | 12 |
| 10-200-52 | Yellow | 100 | Small | 12 |
| 02-539-15 | Blue | 150 | Medium | 12 |
| 02-539-16 | Green | 150 | Medium | 12 |
| 10-200-10 | Red | 150 | Medium | 12 |
| 02-555-109 | Pink | 150 | Medium | 12 |
| 02-539-18 | Yellow | 150 | Medium | 12 |
| 10-200-12 | Blue | 250 | Medium | 12 |
| 10-200-13 | Green | 250 | Medium | 12 |
| 10-200-14 | Red | 250 | Medium | 12 |
| 02-555-110 | Pink | 250 | Medium | 12 |
| 10-200-15 | Yellow | 250 | Medium | 12 |
| 10-200-16 | Blue | 400 | Large | 12 |
| 10-200-17 | Green | 400 | Large | 12 |
| 10-200-18 | Red | 400 | Large | 12 |
| 02-555-111 | Pink | 400 | Large | 12 |
| 10-200-19 | Yellow | 400 | Large | 12 |
| 10-200-20 | Blue | 600 | Large | 6 |
| 10-200-21 | Green | 600 | Large | 6 |
| 10-200-22 | Red | 600 | Large | 6 |
| 02-555-112 | Pink | 600 | Large | 6 |
| 10-200-23 | Yellow | 600 | Large | 6 |
| KIMBLE® ValueWare® Beakers | | | | |
| S00004 | Low Form Griffin | 30 | Small | 48 |
| S00005 | Low Form Griffin | 50 | Small | 48 |
| S00024 | Tall Form | 100 | Small | 48 |
| S00006 | Low Form Griffin | 100 | Small | 48 |
| S00016 | Heavy Duty | 150 | Medium | 48 |
| S00007 | Low Form Griffin | 150 | Medium | 48 |
| S00026 | Tall Form | 200 | Medium | 48 |
| S00017 | Heavy Duty | 250 | Medium | 48 |
| S00028 | Tall Form | 300 | Medium | 48 |
| S00029 | Tall Form | 400 | Large | 36 |
| S00018 | Heavy Duty | 400 | Large | 48 |
| S00009 | Low Form Griffin | 400 | Large | 48 |
| S00031 | Tall Form | 600 | Large | 24 |
| S00019 | Heavy Duty | 600 | Large | 36 |
| S00010 | Low Form Griffin | 600 | Large | 36 |
| S00033 | Tall Form | 1,000 | Large | 18 |
| KIMBLE® KIMAX® Beakers | | | | |
| 02-539-7 | With Spout | 100 | Small | 48 |
| 10-310-5 | No Spout | 100 | Small | 12 |
| 02-539-8 | With Spout | 200 | Medium | 48 |
| 10-310-7 | No Spout | 200 | Medium | 12 |
| 02-555-20A | — | 250 | Medium | 48 |
| 02-539-9 | With Spout | 300 | Medium | 48 |
| 10-310-8 | No Spout | 300 | Medium | 12 |



| Cat. No. | Description | Capacity (mL) | Size Compatibility | Case Qty |
|-------------------------------------------|-------------|---------------|--------------------|----------|
| KIMBLE® KIMAX® Beakers (continued) | | | | |
| 02-539-6 | With Spout | 400 | Large | 36 |
| 02-555-20B | — | 400 | Large | 48 |
| 02-539-11 | With Spout | 600 | Large | 24 |
| 10-310-11 | No Spout | 600 | Large | 6 |
| 02-555-20C | — | 600 | Large | 36 |
| 02-539-12 | With Spout | 1,000 | Large | 18 |
| 10-310-6 | No Spout | 1,000 | Large | 6 |
| KIMBLE® KIMAX® Griffin Beakers | | | | |
| 02-539G | — | 50 | Small | 48 |
| 02-539H | — | 100 | Small | 48 |
| 02-539J | — | 150 | Medium | 48 |
| 02-539K | — | 250 | Medium | 48 |
| 02-539L | — | 400 | Large | 48 |
| 02-539M | — | 600 | Large | 36 |
| 02-539N | — | 800 | Large | 24 |

LOW FORM GRIFFIN BEAKER STARTER PACK

(Consists of one each - 50 mL, 100 mL, 250 mL, 600 mL, 1,000 mL)

| Cat. No. | Size Compatibility | Case Qty |
|----------|-------------------------|----------|
| 02-555-2 | Small, Medium and Large | 1 |

HEAVY DUTY BEAKER STARTER PACK

(Consists of one each - 250 mL, 400 mL, 600 mL, 1,000 mL)

| Cat. No. | Size Compatibility | Case Qty |
|------------|-------------------------|----------|
| 02-555-101 | Small, Medium and Large | 1 |

KIMBLE® COMPLETE SET GLASSWARE STARTER PACK

(Media Bottles, Beakers, Erlenmeyer Flasks, Graduated Cylinders)

| Cat. No. | Size Compatibility | Case Qty |
|------------|-------------------------|----------|
| 12-141-460 | Small, Medium and Large | 1 |

ERLENMEYER FLASKS

| Cat. No. | Description | Capacity (mL) | Size Compatibility | Case Qty |
|---------------------------------------------------|--------------|---------------|--------------------|----------|
| KIMBLE® KIMAX® Erlenmeyer Flasks | | | | |
| 10-091B | Wide Mouth | 250 | Small | 48 |
| 10-091C | Wide Mouth | 500 | Small | 36 |
| 10-310-88 | — | 500 | Small | 1 |
| 10-039K | Narrow Mouth | 1,000 | Small | 24 |
| 10-091D | Wide Mouth | 1,000 | Small | 24 |
| 10-310-88B | — | 1,000 | Small | 1 |
| 10-039M | Narrow Mouth | 2,000 | Small | 8 |
| 10-091E | Wide Mouth | 2,000 | Large | 8 |
| 10-310-88C | — | 2,000 | Small | 1 |
| 10-039P | Narrow Mouth | 4,000 | Small | 1 |
| 10-039Q | Narrow Mouth | 6,000 | Small | 1 |
| KIMBLE® ValueWare® Erlenmeyer Flasks | | | | |
| S00151 | Wide Mouth | 250 | Small | 48 |
| S00152 | Wide Mouth | 500 | Small | 36 |
| S00145 | Narrow Mouth | 500 | Small | 36 |
| S00153 | Wide Mouth | 1,000 | Small | 24 |
| S00146 | Narrow Mouth | 1,000 | Small | 24 |
| S00147 | Narrow Mouth | 2,000 | Small | 8 |
| KIMBLE® KIMAX® Colorware Erlenmeyer Flasks | | | | |
| 10-200-42 | Red | 500 | Small | 6 |
| 10-200-45 | Green | 1,000 | Small | 6 |
| 10-200-32 | Blue | 1,000 | Small | 6 |
| 10-200-46 | Red | 1,000 | Small | 6 |
| 10-200-47 | Yellow | 1,000 | Small | 6 |
| 02-555-118 | Pink | 1,000 | Small | 6 |



ERLENMEYER STARTER PACK

(Consists of one each - 50 mL, 125 mL, 250 mL, 500 mL, 1,000 mL)

| Cat. No. | Size Compatibility | Case Qty |
|-----------|--------------------|----------|
| 02-543-21 | Small | 1 |

GL 45 BOTTLES

| Cat. No. | Description | Capacity (mL) | Size Compatibility | Case Qty |
|------------------------------------------------------|-------------|---------------|--------------------|----------|
| KIMBLE® GL 45 Media Bottles | | | | |
| 10-200-70 | Without Cap | 100 | Small | 10 |
| 02-542A | With Cap | 100 | Small | 10 |
| 02-542B | With Cap | 250 | Small | 10 |
| 02-542C | With Cap | 500 | Small | 10 |
| 02-542D | With Cap | 1,000 | Small | 10 |
| 02-542E | With Cap | 2,000 | Small | 4 |
| 06-421-30 | With Cap | 5,000 | Small | 1 |
| 10-200-76 | Without Cap | 10,000 | Small | 1 |
| 06-421-31 | With Cap | 10,000 | Small | 1 |
| KIMBLE® KIMCOTE® GL 45 Media Bottles | | | | |
| 12-141-320 | — | 100 | Small | 4 |
| 12-141-321 | — | 250 | Small | 4 |
| 12-141-322 | — | 500 | Small | 4 |
| 12-141-323 | — | 1,000 | Small | 4 |
| 12-141-324 | — | 2,000 | Small | 4 |
| 12-141-325 | — | 5,000 | Small | 1 |
| 12-141-326 | — | 10,000 | Small | 1 |
| KIMBLE® RAY-SORB® GL 45 Media Bottles | | | | |
| 02-542-51 | — | 100 | Small | 1 |
| 02-542-52 | — | 250 | Small | 1 |
| 02-542-4C | — | 500 | Small | 1 |
| 02-542-4D | — | 1,000 | Small | 1 |
| 02-542-40 | — | 2,000 | Small | 1 |
| 02-542-41 | — | 5,000 | Small | 1 |
| 02-542-42 | — | 10,000 | Small | 1 |
| KIMBLE® KIMAX® Screw Thread Erlenmeyer Flasks | | | | |
| 10-041-30 | — | 250 | Small | 6 |
| 10-041-31 | — | 1,000 | Small | 6 |
| 10-041-32 | — | 2,000 | Small | 2 |
| WHEATON® LAB 45™ Graduated Media Bottles | | | | |
| 02-940-23 | With Cap | 250 | Small | 12 |
| 02-940-27 | Without Cap | 250 | Small | 12 |
| 06-451-291 | With Cap | 500 | Small | 12 |



| Cat. No. | Description | Capacity (mL) | Size Compatibility | Case Qty |
|-------------------------------------------------------------|-------------|---------------|--------------------|----------|
| WHEATON® LAB 45™ Graduated Media Bottles (continued) | | | | |
| 02-940-32 | Without Cap | 500 | Small | 12 |
| 06-451-292 | With Cap | 1,000 | Small | 12 |
| 02-940-33 | Without Cap | 1,000 | Small | 12 |
| 02-940-24 | With Cap | 2,000 | Small | 6 |
| 02-940-28 | Without Cap | 2,000 | Small | 6 |

KIMBLE® GL 45 MEDIA BOTTLE STARTER PACK

(Pack consists of - 2: 100 mL, 3: 250 mL, 3: 500 mL, 2: 1,000 mL)

| Cat. No. | Size Compatibility | Case Qty |
|------------|--------------------|----------|
| 02-542-355 | Small | 1 |

CYLINDERS

| Cat. No. | Description | Capacity (mL) | Size Compatibility | Case Qty |
|---------------------------------------------------------|---------------------------|---------------|--------------------|----------|
| KIMBLE® KIMAX® TD Economical Graduated Cylinders | | | | |
| 08-549-2F | White Scale, Glass Base | 250 | Small | 18 |
| 08-549-2G | White Scale, Glass Base | 500 | Small | 12 |
| 08-549-2H | White Scale, Glass Base | 1,000 | Medium | 8 |
| 08-549-2J | White Scale, Glass Base | 2,000 | Large | 6 |
| KIMBLE® KimCote® Class A Graduated Cylinders | | | | |
| 12-141-328 | TD | 250 | Small | 1 |
| 12-141-329 | TD | 500 | Small | 1 |
| 12-141-330 | TD | 1,000 | Medium | 1 |
| 12-141-331 | TD | 2,000 | Large | 1 |
| KIMBLE® KIMAX® Class B Graduated Cylinders | | | | |
| 08-549-5G | TD, Single White Scale | 250 | Small | 12 |
| 08-549G | TD, Single Blue Scale | 250 | Small | 12 |
| 08-549-5H | TD, Single White Scale | 500 | Small | 8 |
| 08-549H | TD, Single Blue Scale | 500 | Small | 8 |
| 08-549-5J | TD, Single White Scale | 1,000 | Medium | 4 |
| 08-549J | TD, Single Blue Scale | 1,000 | Medium | 4 |
| 08-549-5K | TD, Single White Scale | 2,000 | Large | 4 |
| 08-549K | TD, Single Blue Scale | 2,000 | Large | 4 |
| KIMBLE® ValueWare® Graduated Cylinders | | | | |
| S00049 | TD, Double White Scale | 250 | Small | 6 |
| S00050 | TD, Double White Scale | 500 | Small | 4 |
| S00051 | TD, Double White Scale | 1,000 | Medium | 4 |
| KIMBLE® ValueWare® Class B Graduated Cylinders | | | | |
| S00041 | TC, White Scale | 250 | Small | 6 |
| S00042 | TC, White Scale | 500 | Small | 4 |
| S00043 | TC, White Scale | 1,000 | Medium | 4 |
| KIMBLE® KIMAX® Class A Graduated Cylinders | | | | |
| 08-548E | TD, Serialized, Certified | 10 to 250 | Small | 4 |
| 08-548-205 | TC | 250 | Small | 4 |
| 08-564A | TD, Reverse Graduations | 10 to 250 | Small | 4 |
| 08-548F | TD, Serialized, Certified | 25 to 500 | Small | 4 |
| 08-548-206 | TC | 500 | Small | 4 |
| 08-567A | TD, Reverse Graduations | 25 to 500 | Small | 4 |
| 08-548G | TD, Serialized, Certified | 50 to 1,000 | Medium | 2 |
| 08-548-207 | TC | 1,000 | Medium | 1 |
| 08-568A | TD, Reverse Graduations | 50 to 1,000 | Medium | 1 |
| 08-548-208 | TC | 2,000 | Large | 1 |




| Cat. No. | Description | Capacity (mL) | Size Compatibility | Case Qty |
|---------------------------------------------------|------------------|---------------|--------------------|----------|
| KIMBLE® KIMAX® Class B Graduated Cylinders | | | | |
| 08-554F | TD, Double Scale | 10 to 250 | Small | 12 |
| 08-566-15C | TC, Single Scale | 250 | Small | 6 |
| 08-554G | TD, Double Scale | 25 to 500 | Small | 8 |
| 08-566-15F | TC, Single Scale | 500 | Small | 4 |
| 08-554H | TD, Double Scale | 50 to 1,000 | Medium | 4 |
| 08-566-15D | TC, Single Scale | 1,000 | Medium | 4 |
| 08-554J | TD, Double Scale | 100 to 2,000 | Large | 4 |
| 08-566-15H | TC, Single Scale | 2,000 | Large | 2 |

RESERVOIR BOTTLES

| Cat. No. | Description | Capacity (mL) | Size Compatibility | Case Qty |
|---------------------------------------------------|----------------------|---------------|--------------------|----------|
| KIMBLE® ULTRA-WARE® Reservoir Bottles | | | | |
| K959902-0502 | Ungraduated | 500 | Small | 1 |
| K953902-0502 | Graduated | 500 | Small | 1 |
| K953901-2002 | Ungraduated | 1,000 | Small | 1 |
| K953902-2002 | Ungraduated | 2,000 | Small | 1 |
| K953902-5002 | Ungraduated | 5,000 | Small | 1 |
| K953902-1003 | Ungraduated | 10,000 | Small | 1 |
| K953901-1003 | Graduated | 10,000 | Small | 1 |
| 22-302625 | Ungraduated | 20,000 | Small | 1 |
| KIMBLE® ULTRA-WARE® HPLC Reservoir Bottles | | | | |
| K953901-0252 | Ungraduated | 250 | Small | 1 |
| K953901-0502 | Ungraduated | 500 | Small | 1 |
| K953901-1002 | Ungraduated | 1,000 | Small | 1 |
| K953901-2002 | Ungraduated | 2,000 | Small | 1 |
| K953901-5002 | Ungraduated | 5,000 | Small | 1 |
| K953910-5002 | Graduated, Side Neck | 5,000 | Small | 1 |
| K953901-1003 | Ungraduated | 10,000 | Small | 1 |
| K953901-2003 | Ungraduated | 20,000 | Small | 1 |





KIMBLE® is a registered trademark of DWK Life Sciences LLC or its subsidiaries.
FBRO_023 12/2019

Contact us today:

In the United States

Order online: fishersci.com
Fax an order: 1-800-926-1166
Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca
Fax an order: 1-800-463-2996
Call customer service: 1-800-234-7437

