



Fisherbrand[®]
QUALITY. RELIABILITY. VALUE.

Fisherbrand Polymer Carts

Carts designed for a wide variety of applications.

Specifications for BC Series and Deep Ledge

Polymer carts are available in 2 or 3 shelf models with either a 1/2" raised ship's edge or 2 3/4" Deep Ledge design and four durable chrome-plated posts. The high density polyethylene shelves offer superior corrosion and chemical resistance. Carts roll effortlessly on four (4) 4" diameter swivel casters. Cart weight capacity is 150 lbs. (68kg.) per shelf, 400 lbs. (181 kg.) per cart.

BC Series Carts

Cart is available in 2-shelf and 3-shelf.

Catalog Number	Weight/Length/Height (in)	Description	Approx. Pkd. Wt. (lbs.)
13451181	18 x 29 1/2 x 33 1/4	3-shelf unit	33
13451182	21 1/2 x 33 3/4 x 33 1/4	2-shelf unit	34 1/2
13451183	21 1/2 x 33 3/4 x 33 1/4	3-shelf unit	43

Deep Ledge Series Carts

Cart is available in 2-shelf and 3-shelf.

Catalog Number	Weight/Length/Height (in)	Description	Approx. Pkd. Wt. (lbs.)
13451184	21 1/2 x 32 3/4 x 41	3-shelf unit	37 1/2
13451185	27 x 38 3/4 x 41	2-shelf unit	35
13451186	27 x 38 3/4 x 41	3-shelf unit	46 1/2

- **Easy, No-Tool Assembly:** Ships knocked down.
- **Ultimate Cleanability:** Smooth surfaces are designed to resist staining and wipe clean easily.
- **Corrosion-Proof Shelves:** All polymer shelf will never corrode.
- **Contains Spills:** BC Series shelves feature 1/2" raised ship's edge design; Deep Ledge Series shelves feature a deeper 2 3/4" ledge to contain product spills.
- **Attractive Aesthetics:** allow for use in a variety of applications.
- **Sturdy Construction:** Each shelf can hold up to 150 lbs.; each unit can hold up to 400 lbs.
- **Easy-to-Adjust:** Optional center shelf adjusts on 1" increments.
- **Ergonomic Handle** designed for easy use.

Fisherbrand is a registered trademark of Fisher Scientific.



3-Shelf Deep Ledge Cart



2-Shelf Deep Ledge Cart



2-Shelf BC Cart



3-Shelf BC Cart



Fisher Scientific

For customer service, call 1-800-766-7000.
To fax an order, use 1-800-926-1166.
To order online: www.fishersci.com