



30+
YEARS
OF GLOVE EXPERTISE

A PERFECT PAIR

CHOOSING THE RIGHT GLOVE FOR YOU

Find it at fishersci.com



WHY CHOOSE MEDICOM GLOVES



Responsible Supply Chain

We are committed to maintaining the highest ethical standards in our supply chain. All our suppliers must sign a code of conduct and pass regular independent audits.



Innovative Glove Design and Development

Our passionate product development teams drive continuous innovation to design high-quality gloves that not only help eliminate pin holes and suit a variety of uses, but also ensure wearer comfort.



Robust Manufacturing Processes and Quality Controls

Our manufacturing processes ensure consistent glove quality and fit. We meet or exceed performance standards for protection and durability.



Uncompromising Material Quality

Our raw materials are carefully selected and tightly controlled to prevent impurities in finished gloves and all steps throughout production are monitored to ensure minimal bioburden.



30+ Years of Glove Expertise

Hand protection remains a key focus within our infection prevention and control portfolio



CHOOSING THE RIGHT GLOVES

Looking for the perfect match? Here's what you need to know when choosing gloves.

Which material best fits your needs?

LATEX

- Highest level of comfort, flexibility and tactile sensitivity
- Optimal fit
- Higher degree of protection against viruses and bacteria

NITRILE

- Optimal latex-free option
- Puncture resistant
- Comfortable fit
- Strong protection against chemicals

VINYL

- Best value latex-free gloves option
- Good fit

POLYCHLOROPRENE

- Comfort of latex
- Strength of nitrile
- Optimal fit

FIT AND FEEL

Texture

Textured fingertips ensure ensure a secure grip, even when wet. Fully textured gloves provide ultimate grip.

Inner treatment

Powder-free gloves cause less irritation and are cleaner to use because they do not leave residue on clothing or surfaces. They are also safer for use in surgeries and other invasive procedures.

Colour

Our gloves are available in a wide range of colours to meet a variety of needs and preferences.

Cuffs

Beaded cuffs provide better coverage and protection for high-risk tasks. They also increase glove strength when donning.



QUICK TIPS

DO YOU WANT A BETTER GRIP?

Textured fingertips provide a more secure grasp, especially when gloves are wet. Many professionals appreciate the added control. For optimal grip, choose one of our fully textured gloves.

THE IMPORTANCE OF FIT

Gloves that are too small can contribute to hand fatigue and cause injuries such as carpal tunnel syndrome. They are also less resistant to punctures and tears. Gloves that are too large decrease tactile sensitivity, reduce manual dexterity and increase the risk of cross-contamination.

SPECIALITY CONSIDERATIONS



CHEMO SAFE

Gloves that are tested and approved for use with chemotherapy drugs



DERMASENSE TECHNOLOGY

Helps reduce the risk of skin reactions that could be caused by chemical accelerators commonly found in nitrile gloves



FENTANYL SAFE

Gloves that are tested and approved for use with fentanyl



CHEMICAL RESISTANT

Gloves that are tested and approved for use with certain chemicals

QUICK TIPS

CONSIDER THE INTENDED USE OF YOUR GLOVES

Surgical procedures require gloves that fit closer to the hand for greater tactile sensitivity. For patient examinations, blood tests and patient care, medical examination gloves are preferred by most healthcare professionals. For handling drugs, chemical components or food, choose gloves that comply with industry standards for these tasks.

DO YOU HAVE A LATEX ALLERGY/SENSITIVITY, CONTACT DERMATITIS OR OTHER SKIN CONDITION?

Vinyl, nitrile or polychloroprene gloves are good options. Powder-free and gloves without chemical accelerators also have a lower risk of causing allergic reactions or other skin reactions.

FIND YOUR MATCH

What criteria are most important to you? Use the chart to find the glove that suits you best.



Order your samples

Product Name	Material	Thickness	Specialty Considerations	Grip	Skin-friendly Considerations	Color
SafeBasics® True Fit® Thin	NITRILE	Thin	Chemo Tested / Fentanyl Safe	Textured Fingertips		Blue
SafeBasics® True Fit® Thin	NITRILE	Thin	Chemo Tested	Textured Fingertips		Pink
AssureTouch™ Basics	NITRILE	Thin	Chemo Tested	Textured Fingertips		Blue
AssureTouch™ Agile™	NITRILE	Medium	Chemo Tested	Textured Fingertips		Blue
AssureTouch™ Power™	NITRILE	Medium	Chemo Tested / Fentanyl Safe	Textured Fingertips		Blue
AssureTouch™ Shield™	VINYL	Medium				Clear
AssureTouch™ Ultra Tough™	NITRILE	Thick		Fully Textured		Black
Ocean Pacific® NeoNatural™	POLY CHLORO	Medium	Chemo Tested / Chemical Resistant	Textured Fingertips		Emerald

Medical Grade Gloves

All our gloves are fully tested and comply with ASTM standards.

GET TO KNOW YOUR GLOVES

You have found the right glove! Now, get to know everything about it.



NITRILE SafeBasics® True Fit™ Thin

Thinner chemo-tested glove with the unrivaled strength and durability of nitrile for reliable protection and exceptional fit. Textured fingertips ensure a secure grip. Item 1185 is also fentanyl safe.

Item		1185	1186
Color		Blue ●	Pink ●
Palm Thickness/Weight		3 mil, 3.5 g	2.5 mil, 3.2 g
Size	XS	17101055	17101060
	S	17101056	17101061
	M	17101057	17101062
	L	17101058	17101063
	XL	17101059	
Packaging		300 gloves/box	



NITRILE AssureTouch™ Basics

This affordable chemo-tested glove provides healthcare professionals with the hand protection they need to prevent exposure and contamination when handling chemotherapy drugs. High-quality nitrile won't compromise manual dexterity or tactile sensitivity. Textured fingertips ensure a secure grip, even when wet.

Item		9985
Color		Blue ●
Palm Thickness/Weight		2.7 mil, 3.5 g
Size	XS	17101088
	S	17101089
	M	17101090
	L	17101091
	XL	17101092
Packaging		100 gloves/box



NITRILE AssureTouch™ Agile™

This chemo-tested glove offers optimal puncture, tear and chemical resistance without compromising dexterity or tactile sensitivity. Textured fingertips ensure a secure grip, even when wet.

Item		9998
Color		Blue ●
Palm Thickness/Weight		4.5 mils, 6.2 g
Size	XS	17101083
	S	17101084
	M	17101085
	L	17101086
	XL	17101087
Packaging		100 gloves/box



POLY CHLORO Ocean Pacific® NeoNatural™

Polychloroprene combines the softness of latex with the strength of nitrile. Textured fingertips, latex tested and chemical resistant.

Item		1228
Color		Emerald ●
Palm Thickness/Weight		4 mil, 6 g
Size	XS	17101050
	S	17101051
	M	17101052
	L	17101053
	XL	17101054
Packaging		100 gloves/box



NITRILE AssureTouch™ Power™

This chemo-tested and fentanyl-safe glove provides strong puncture, tear and chemical resistance without compromising dexterity or tactile sensitivity. Extended cuffs provide optimal coverage. Textured fingertips for a secure grip, even when wet.

Item		9999
Color		Blue ●
Palm Thickness/Weight		4 mil, 6.2 g
Size	XS	17101073
	S	17101074
	M	17101075
	L	17101076
	XL	17101077
Packaging		100 gloves/box



VINYL AssureTouch™ Shield™

Versatile vinyl glove with a smooth texture for enhanced tactile sensitivity.

Item		9994
Color		Clear ●
Palm Thickness/Weight		3.2 mil, 5.3 g
Size	S	17000809
	M	17000810
	L	17000811
	XL	17000812
Packaging		100 gloves/box



Medicom[®]
Pride in Protection[™]

For more information, please visit medicom.com or contact us at:
cservice@medicom.ca / 1-800-361-2862

AMD Medicom Inc. 2555 Chemin de l'Aviation, Montreal, Quebec
H9P 2Z2 Canada

©2023 AMD Medicom Inc. All rights reserved.

Rev. 01

Distributed by Fisher Scientific. Contact us today:

In the United States

Order online: fishersci.com

Call customer service: 1-800-766-7000

 **fisher scientific**
part of **Thermo Fisher Scientific**

© 2023 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at fishersci.com/trademarks.

BN225340108-US