



AssureMask[®] Balance[™]

Procedure Earloop Face Masks



The perfect balance between comfort and protection

AssureMask Balance earloop face masks provide healthcare professionals and other frontline workers with the ideal balance between the protection they need and the breathability and comfort they want. Made from high-quality materials, they protect without the gaping, discomfort or irritation often caused by masks made from lower-quality materials.

Available colors:    



Find it at fishersci.com and fishersci.ca



PROTECTION. COMFORT. BREATHABILITY.



Shingle pleat design
prevents fluid pooling and
cross-contamination



Made from **high-quality materials** that provide
reliable protection without
compromising breathability



Available in
ASTM levels 1, 2 and 3



Adjustable nosepiece
ensures a comfortable,
customized fit



Ultrasonically sealed,
soft, stretchy, earloops
free of natural rubber
latex prevent pulling and
irritation

ASTM
LEVEL 1

19910757	White
19910758	Blue*
19910759	Pink
19910760	Yellow

50 masks/box
10 boxes/case

*Canada Catalog No.:296368455

ASTM
LEVEL 2

17000545	White
19910740	Blue
17000546	Pink
17000271	Yellow

50 masks/box
10 boxes/case

ASTM
LEVEL 3

17000547	White
17000544	Blue
17000548	Pink
17000549	Yellow

50 masks/box
10 boxes/case

Distributed by Fisher Scientific. Contact us today:

In the United States

Order online: fishersci.com

Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca

Call customer service: 1-800-234-7437

 **fisher scientific**
part of Thermo Fisher Scientific