



AssureTouch™ Ultra Tough

Nitrile Medical Examination Gloves



EXTRA PROTECTION AND GRIP when you need it

Made from industrial quality nitrile to provide strong tear resistance without compromising comfort or fit. **Fully textured** to ensure a secure grip, even when wet.



Find it at fishersci.com and fishersci.ca





KEY FEATURES



Fully textured



Tear resistant



Comfortable fit

BENEFITS OF NITRILE GLOVES

Comfortable fit

Nitrile gloves conform to hands, providing a snug, comfortable fit and reducing hand fatigue.

Optimal stretch

Nitrile offers optimal flexibility, making nitrile gloves ideal for procedures that require manual dexterity.

Latex free

Nitrile gloves are not made with natural rubber latex, so they are safe for professionals and patients with a latex allergy or sensitivity.



17101093 Small	17101096 Extra Large
17101094 Medium	17101097 2X Large
17101095 Large	

5 mil, 6.5 g
100 gloves/box
10 boxes/case

Distributed by Fisher Scientific. Contact us today:

In the United States

Order online: fishersci.com
Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca
Call customer service: 1-800-234-7437

 **fisher scientific**
part of Thermo Fisher Scientific