

AssureWear™

Over-the-head Isolation Gown



**COVERAGE,
COMFORT AND
CONVENIENCE**
when it matters most

AAMI PB70 Level 2 protection in an isolation gown that's easy to don and doff.

- Ultrasonically sealed, over-the-head neck closure makes it easy to put on and take off when every second counts.
- Full back closure provides complete coverage for optimal protection.
- Flat belt on the left-hand side creates a personalized fit.
- Elastic thumb loop cuffs prevent sleeves from riding up.

The AAMI-rated protection you need in a gown that's as easy to don and doff as it is to wear.

Features	Benefits
Ultrasonically sealed over-the-head neck closure	Allows safe and easy donning and doffing
AAMI PB70 level 2 rating	Offers effective barrier performance for low fluid level procedures
Full back closure	Provides complete coverage for optimal protection
Elastic, thumb loop cuffs	Prevent sleeves from riding up
Non-woven SMMS fabric	Keeps wearers cool and comfortable
Flat belt on left hand side	Ensures a personalized fit

17000575 Large, Yellow

17000576 Extra Large, Yellow

10 gowns per bag; 10 bags per case; 100 gowns per case



Trust

With over 30 years of infection control expertise, Medicom is one of the world's most trusted manufacturers and distributors of personal protective equipment (PPE) and infection control solutions.



Quality

Medicom gowns are manufactured on state-of-the-art machinery in our own factories. Our innovative designs and premium materials ensure reliable protection and exceptional comfort and breathability.



Protection

In addition to gowns, Medicom offers a broad range of medical masks and other personal protective equipment products. **Discover our complete portfolio at [medicom.com](https://www.medicom.com).**

Distributed by Fisher Scientific. Contact us today:

In the United States

Order online: fishersci.com
Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca
Call customer service: 1-800-234-7437

