

TIRED OF INTERRUPTIONS?

A Better Solution Has Arrived

WHEATON® μ L MicroLiter® BlueMAG™ Screw Caps

- ✓ Safer, more robust magnetic performance
Guards against dropped or misaligned vials
- ✓ Better grip reduces closure errors and speeds up work
Run automation at maximum speed!
- ✓ Fits 9mm, 16mm and 18mm vials for HPLC/LCMS and GC/GCMS use
Multiple septa options!
- ✓ Universal compatibility with lab automation equipment
Prep work stations, Autosamplers, Magnetic X-Y-Z robots!



9mm BlueMAG™ Polypropylene Screw Caps — Ideal for LC-MS

Cat. No.	Description	Liner	Qty / Case
15-112-109	Screw Cap with Steel Overlay	White Silicone /Red PTFE; 0.040"	100
15-112-110	Screw Cap with Steel Overlay	White Silicone /Red PTFE; 0.040", single slit	100
15-112-111	Screw Cap with Steel Overlay	White Silicone /Red PTFE/White Silicone; 0.040"	100



18mm BlueMAG™ Polypropylene Screw Caps — Ideal for GC and GC-MS applications

NOTE: 18mm BlueMAG cap fits both 16mm and 18mm headspace vials

Cat. No.	Description	Liner	Qty / Case
15-112-102	Screw Cap with Steel Overlay	White Silicone/Blue PTFE Microcenter SPME	100
15-112-103	Screw Cap with Steel Overlay	White Silicone/Natural PTFE; 0.125"	100
15-112-104	Screw Cap with Steel Overlay	Red Butyl/Natural PTFE; 0.125"	100
15-112-105	Screw Cap with Steel Overlay	Trans Blue Silicone /White PTFE; 0.070"	100
15-112-106	Screw Cap with Steel Overlay	White Silicone /Blue PTFE; 0.060"	100
15-112-107	Screw Cap with Steel Overlay	White Silicone /Red PTFE; 0.075"	100
15-112-108	Screw Cap with Steel Overlay	Red Rubber /Natural PTFE; 0.065"	100



Contact us today:

In the United States

Order online: fishersci.com
Fax an order: 1-800-926-1166
Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca
Fax an order: 1-800-463-2996
Call customer service: 1-800-234-7437



fisher scientific
part of Thermo Fisher Scientific

BN20191362