




DWK LIFE SCIENCES
DURAN WHEATON KIMBLE

WHEATON® CENTRIFUGE TUBES

Ideal for standard laboratory sample management
centrifugation needs



Find it at [fishersci.com](https://www.fishersci.com) and [fishersci.ca](https://www.fishersci.ca)

 **fisher scientific**
part of Thermo Fisher Scientific

WHEATON® POLYPROPYLENE CENTRIFUGE TUBES

IDEAL FOR STANDARD LABORATORY CENTRIFUGATION NEEDS



WHEATON® Centrifuge Tubes

Manufactured from ultra-clear virgin polypropylene, the WHEATON® Centrifuge Tube range meets the needs of most laboratory users for routine centrifugation applications. Available in sizes 15mL and 50mL, with a leak-free cap and large writing area, they are ideal for a wide range of sample handling applications.

- Manufactured from virgin ultra-clear polypropylene for good sample visibility
- Flat cap design:
 - Offers excellent leak-free performance
 - Easy to write on and mark samples
- White-solvent-resistant printed graduations with large marking strip for labelling
- RCF values up to 9500xg
- Non-pyrogenic to < 0.25EU/mL
- RNase and DNase free
- Bulk packed in a peel open-bag
- Supplied ETO sterile
- Easy to handle with gloved hands

PEEL BAG PACK



Cat. No.	Description	Max RCF (xg)	Graduation Range (ml)	Graduation Intervals (ml)	Shelf Pack Qty	Case Qty
15112202	15ml Centrifuge Tube, Polypropylene, Bagged, Sterile	9000	2 to 14.5	0.5	50	500
15112201	50ml Centrifuge Tube, Polypropylene, Bagged, Sterile	9500	5 to 47.5	2.5	25	500
15112203	50ml Centrifuge Tube, Polypropylene, Self-Standing, Bagged, Sterile	7000	5 to 47.5	2.5	25	500

Contact us today:

In the United States

Order online: fishersci.com

Fax an order: 1-800-926-1166

Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca

Fax an order: 1-800-463-2996

Call customer service: 1-800-234-7437

 **fisher scientific**
part of **Thermo Fisher Scientific**