



DWK LIFE SCIENCES
DURAN WHEATON KIMBLE

WHEATON® PLASTIC PASTEUR PIPETTES

Ideal for the accurate volume transfer of liquid
and bacteriological samples.



Find it at [fishersci.com](https://www.fishersci.com) and [fishersci.ca](https://www.fishersci.ca)

 **fisher scientific**
part of Thermo Fisher Scientific

WHEATON® SINGLE-USE POLYETHYLENE PASTEUR PIPETTES

IDEAL FOR GENERAL LIQUID TRANSFER APPLICATIONS AND
EXTRACTING SAMPLES FOR ANALYSIS



WHEATON® Plastic Pasteur Pipettes

WHEATON® polyethylene pasteur pipettes are an everyday laboratory essential. Available in a range of sizes from 1 mL graduated to a jumbo-sized to meet user applications for use in areas such as blood banking, bacteriology, general life science and industry. The integral bulb design and single-use feature eliminate the need for plugging or autoclaving, helping to reduce the potential of cross-contamination during everyday sample handling.

- Manufactured from soft, flexible, inert non-toxic polyethylene
- Ideal for safe and easy transfer of liquid samples
- Clear graduations on the 1 mL and 3 mL versions support ease of volume control
- Suitable for use as a freeze vial when the bulb is fully heat sealed
- Jumbo pipettes are ideal for larger volume transfer or accessing liquids in taller vessels
- Easy to handle with gloved hands
- Non-toxic material can be safely incinerated for disposal

Cat. No.	Description	Sterility	Graduated	Graduation Intervals (mL)	Overall Length (cm)	Working Capacity (mL)	Bulb Draw (mL)	Drop Size (uL)	Shelf Pack Qty	Case Qty
15112204	1 mL Graduated Pasteur Pipette, Bulk	NS	up to 1 mL	0.25	15	5	3.5	34.5	500 x 6 bags	3000
15112205	3 mL Graduated Pasteur Pipette, Bulk	NS	up to 3 mL	0.5	16	7	3.5	47.6	500 x 6 bags	3000
15112206	Jumbo Pasteur Pipette, Bulk	NS	N/A	N/A	30	23	10	22	100 x 6 bags	600

Contact us today:

In the United States

Order online: fishersci.com

Fax an order: 1-800-926-1166

Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca

Fax an order: 1-800-463-2996

Call customer service: 1-800-234-7437

