



Nitrile Exam Gloves, Extended Cuff

- UL Approved for Emergency Medical Use
- Chemotherapy approved
- Thickness: 6.0 mil (palm), 7.0 mil (finger)
- Finish: Textured
- Color: Teal
- Length: 12" Extended Cuff
- Non-sterile/Ambidextrous
- 100% Latex-free
- 50 Gloves/Box, 10 Boxes/Case

Size	Cat. No.
Small	19-041-170A
Medium	19-041-170B
Large	19-041-170C
X-Large	19-041-170D



Comfort Nitrile Exam Gloves

- Enhanced Comfort and Sensitivity
- Thickness: 3.0 mil (palm), 3.5 mil (finger)
- Finish: Textured Fingers
- Color: Periwinkle Blue
- Length: 9.5" Cuff
- Non-sterile/Ambidextrous
- 100% Latex-free
- 200 Gloves/Box, 10 Boxes/Case
Size XL contains 180 Gloves/Box

Size	Cat. No.
X-Small	19-041-171A
Small	19-041-171B
Medium	19-041-171C
Large	19-041-171D
X-Large	19-041-171E



Comfort Nitrile Exam Gloves with Aloe

- Soothing Aloe Surrounds the Hands
- Thickness: 3.0 mil (palm), 4.0 mil (finger)
- Finish: Textured
- Color: Blue Green
- Length: 9.5" Cuff
- Non-sterile/Ambidextrous
- 100% Latex-free
- 100 Gloves/Box, 10 Boxes/Case

Size	Cat. No.
X-Small	19-041-172A
Small	19-041-172B
Medium	19-041-172C
Large	19-041-172D
X-Large	19-041-172E



Nitrile Exam Gloves Powder-Free

- Thickness: 3.0 mil (palm), 4.3 mil (finger)
- Finish: Textured
- Color: Blue
- Length: 9.5" Cuff
- Non-sterile/Ambidextrous
- 100% Latex-free
- 100 Gloves/Box, 10 Boxes/Case

Size	Cat. No.
X-Small	19-130-1597A
Small	19-130-1597B
Medium	19-130-1597C
Large	19-130-1597D
X-Large	19-130-1597E
XX-Large	19-130-1597F



**Latex Exam Gloves
Powder-Free**

- Thickness: 5.0 mil (palm), 6.3 mil (finger)
- Finish: Textured
- Color: Natural
- Length: 9.5"
- Non-sterile/Ambidextrous
- Low-protein 510K rating
- 100 Gloves/Box, 10 Boxes/Case

Size Cat. No.

Small	11-394-5A
Medium	11-394-5B
Large	11-394-5C
X-Large	11-394-5D



**Latex Exam Gloves
Powdered**

- Thickness: 5.0 mil (palm), 6.0 mil (finger)
- Finish: Smooth
- Color: Natural
- Length: 9.5"
- Non-sterile/Ambidextrous
- Low-protein 510K rating
- 100 Gloves/Box, 10 Boxes/Case

Size Cat. No.

X-Small	11-394-4AA
Small	11-394-4A
Medium	11-394-4B
Large	11-394-4C
X-Large	11-394-4D



**Aloe Nitrile Exam Gloves
Powder-Free**

- Chemotherapy approved
- Soothing aloe coating
- Thickness: 5.5 mil (palm), 6.0 mil (finger)
- Finish: Textured
- Color: Green
- Length: 9.5"
- Non-sterile/Ambidextrous
- 100% latex-free
- 100 Gloves/Box, 10 Boxes/Case

Size Cat. No.

Small	19-050-550A
Medium	19-050-550B
Large	19-050-550C
X-Large	19-050-550D



**Aloe Latex Exam Gloves
Powder-Free**

- Thickness: 6.0 mil (palm), 6.5 mil (finger)
- Finish: Textured
- Color: Natural
- Length: 9.5"
- Non-sterile/Ambidextrous
- Low-protein 510K rating
- 100 Gloves/Box, 10 Boxes/Case

Size Cat. No.

Small	19-050-548A
Medium	19-050-548B
Large	19-050-548C
X-Large	19-050-548D

WARNING: This product contains natural rubber latex which may cause allergic reactions. A small percentage of the population may have a heightened sensitivity to natural rubber latex, and prolonged use may cause allergic reactions in such persons. In the event you have, or suspect you may have, such hypersensitivity, you are advised to seek non-latex alternatives. If during use rashes or other signs of discomfort occur, discontinue use immediately and consult your physician. Safe use of this glove by or on latex-sensitized individuals has not been established.

In the United States:

For customer service, call 1-800-766-7000
To fax an order, use 1-800-926-1166
To order online: www.fishersci.com