



Fisherbrand[®]
QUALITY. RELIABILITY. VALUE.

Fisherbrand Carts

Designed with your lab in mind.

Overview

Fisherbrand Laboratory Carts fit in every lab. These versatile carts have endless uses from transporting equipment and supplies to providing auxiliary bench space.

Catalog No.	L/W/H (in)	Description	Wt (lbs)
15252400	34.8 x 19 x 36.4	Stainless Steel Cart. Features two SST pans that are 2" deep on three sides with one open side for easy transfer of equipment from countertop to cart	56
15252401	34.8 x 19 x 36.4	2-Basket Glassware Cart. Used to carry glassware and small instruments in 2 vinyl-coated wire baskets. Includes removable plastic drip pan to catch spills	43
15252402	34.8 x 19 x 36.4	4-Basket Glassware Cart. Used to carry glassware and small instruments in 4 vinyl-coated wire baskets. Includes removable plastic drip pan to catch spills	43
15252403	34.8 x 19 x 36.4	Pan Cart. Used for transport of wet or soiled labware. Features 2" deep powder coated steel, removable watertight pans lined with gray rubber mats to prevent shifting	58
15252404	34.8 x 19 x 36.4	Portable Table. Features two open shelves that provide easy transfer of equipment from countertop to cart	44
15252405	34.8 x 19 x 36.4	Chemical Cart. Features two open shelves with rigid PVC extrusions on two sides to prevent containers from sliding	46

Fisherbrand is a registered trademark of Fisher Scientific.

All Carts Feature:

- 4" diameter casters, two with toe locks
- Ready-to-roll, full factory assembly
- Ergonomic, extruded rubber handle grips
- Durable Fisherbrand blue, powder-coated 18 gauge welded tubular steel frame
- ID tag
- 5 year warranty
- Supports loads up to 400 lbs.



15252401



15252402



15252400



15252403



15252404



15252405



Fisher Scientific

For customer service, call 1-800-766-7000.
To fax an order, use 1-800-926-1166.
To order online: www.fishersci.com

© 2013 Thermo Fisher Scientific Inc. All rights reserved.

BN0103133