



Ocean Pacific® Ninja™

Nitrile Medical Examination Gloves



EXCEPTIONAL PROTECTION

when you need it most

These thicker, **fully textured** black nitrile gloves provide optimal grip, puncture resistance and protection against bloodborne pathogens, chemical lubricants and sharp instruments. **Chemo tested and fentanyl safe.**



Find it at fishersci.com and fishersci.ca

 **fisher scientific**
part of Thermo Fisher Scientific



**"These thicker, fully textured black gloves provide optimal grip and protection against sharp instruments."
- Medicom R&D Team**

KEY FEATURES



Fully textured



Fentanyl safe



Chemo tested

BENEFITS OF NITRILE GLOVES

Comfortable fit

Nitrile gloves conform to hands, providing a snug, comfortable fit to reduce hand fatigue.

Optimal stretch

Nitrile gloves offer optimal flexibility, making them ideal for procedures that require manual dexterity and tactile sensitivity.

Latex free

Nitrile gloves are not made with natural rubber latex, so they are safe for professionals and patients with a latex allergy or sensitivity.



17101069 Small

17101071 Large

17101070 Medium

17101072 Extra Large

Palm thickness/weight: 4.8 mil

Weight: 5.3 g

200 gloves/box

Distributed by Fisher Scientific. Contact us today:

In the United States

Order online: fishersci.com

Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca

Call customer service: 1-800-234-7437

 **fisher scientific**
part of Thermo Fisher Scientific