



# SafeMask<sup>®</sup> Master Series

Procedure Earloop Face Masks

**BURST**  
of colour



Engineered with  
**Simply Soft Technology™**

Find it at [fishersci.com](https://fishersci.com) and [fishersci.ca](https://fishersci.ca)



# SafeMask® Master Series

Procedure Earloop Face Masks

## Colour Your World

- A beautiful line of masks designed in a variety of unique and vibrant colours so that your mask is as personalized as the dental care you provide.
- Master Series ASTM Masks are engineered with Simply Soft Technology™ for a smooth, luxurious and comfortable fit, a perfect blend of breathability and softness.

ASTM LEVEL 1		ASTM LEVEL 2		ASTM LEVEL 3	
<b>19910661</b> Augusta Sky Royal Blue	<b>19910663</b> Lush Lawn Vibrant Green	<b>19910667</b> Ocean Surf Aquamarine	<b>19910666</b> Azalea Festival Bright Fuchsia	<b>19910668</b> Ocean Surf Aquamarine	<b>19910669</b> Azalea Festival Bright Fuchsia
<b>19910662</b> Azalea Festival Bright Fuchsia	<b>19910665</b> Southern Bellflower Radiant Orchid	<b>19910666</b> Azalea Festival Bright Fuchsia		<b>19910670</b> Southern Bellflower Radiant Orchid	

### Performance factors

	ASTM LEVEL 1	ASTM LEVEL 2	ASTM LEVEL 3
Bacterial filtration efficiency (BFE)	≥ 95%	≥ 98%	≥ 98%
Sub-micron particulate efficiency at 0.1 µm (PFE)	≥ 95%	≥ 98%	≥ 98%
Resistant to penetration by synthetic blood	@80 mm Hg	@120 mm Hg	@160 mm Hg
Differential Pressure (ΔP)	< 4.0 mm H <sub>2</sub> O/cm <sup>2</sup>	< 5.0 mm H <sub>2</sub> O/cm <sup>2</sup>	< 5.0 mm H <sub>2</sub> O/cm <sup>2</sup>
Flame	Class 1	Class 1	Class 1

Suggested use  
Packaging configuration

Daily favourite for everyday procedures  
50 masks per box, 10 boxes per case



Made in  
 USA

Distributed by Fisher Scientific. Contact us today:

**In the United States**

Order online: fishersci.com  
Call customer service: 1-800-766-7000

**In Canada**

Order online: fishersci.ca  
Call customer service: 1-800-234-7437

**fisher scientific**  
part of Thermo Fisher Scientific