



**PREFERRED
VENDOR PRICING**

GILSON ON-SITE PIPETTE SERVICE

Never question the integrity of your samples!

Our flexible, tiered service level offerings, along with the option for ISO 17025 accredited service, are tailored to fit your lab's needs, and designed to minimize unplanned, costly repairs while increasing the lifespan of your pipettes.

WE SERVICE ALL BRANDS

📍 **PLACE ·**
DATES ·

CONTACT US TO SCHEDULE YOUR APPOINTMENT



Find it at fishersci.com and fishersci.ca

 **fisher scientific**
part of Thermo Fisher Scientific

HOW TO ORDER SERVICE GUIDE

WE SERVICE ALL BRANDS OF PIPETTES

GILSON PERFORMED SERVICE ON-SITE

- Contact your Fisher Rep to sign up for the service event.
- On the day of the event, take your pipettes to the Gilson Service Technician. Fill out the COF authorizing Gilson to work on your pipettes. (provided by Gilson and Fisher).
- You will be notified when your pipettes are ready to be picked up.
- You will be provided with quantity and part numbers to place your order with Fisher.
- Place your order with Fisher prior to picking up your pipettes. Once your order is placed you will be able to pick up your pipettes.

WHAT'S INCLUDED

LITE (Perfect For Academia)
DEPENDABLE PREVENTATIVE MAINTENANCE

PREMIUM
OUR MOST COMPREHENSIVE SERVICE OFFERING

Visual & Functional Inspection	◆	◆
Internal & External Cleaning	◆	◆
Leak Test	◆	◆
Gravimetric Check	◆	◆
Preventative Maintenance	◆	◆
Re-grease As Needed	◆	◆
Pipette Label	◆	◆
Mail-In or Onsite Service Available	◆	◆
Calibration & Detailed Certificate		◆
ISO 17025 Calibration Available		◆
Lifetime Warranty		◆
Free Repairs		◆

* Includes all replacement parts at the discretion of the service technician. Parts need to be associated with a reasonable cost. See Gilson Terms of Pipette Service on gilson.com.

PART NUMBER	PRICE (PER PIPETTE)	DESCRIPTION
SVC101199G	\$36	PREMIUM Single Channel Pipette Service (Repairs Included)
SVC101200G	\$44	PREMIUM Multichannel Service (Repairs Included)
SVC101101G	\$23	LITE Single Channel Pipette Service with Balance Check
SVC101102G	\$33	LITE Multichannel Pipette Service with Balance Check

CONTACT US TO SCHEDULE YOUR APPOINTMENT

Distributed by Fisher Scientific. Contact us today:

In the United States

Order online: fishersci.com

Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca

Call customer service: 1-800-234-7437



24GIL4359

